



∴ **HEINE GAMMA® Series**


Shockproof and Latex-free Sphygmomanometers
and Stethoscopes



∴ Best in Class.

Those who expect the highest diagnostic precision, choose the best quality.


- ∴ Shockproof.
- ∴ Double-walled housing.
- ∴ Three-point manometer mounting.
- ∴ Exceeds all relevant standards.



Double-walled housing with three-point manometer mounting


GAMMA G Series is built tough, incorporating an intricate shock-absorbing system that protects the manometer and produces instruments of unsurpassed durability. The HEINE GAMMA G Series is tested well above the European (EN1060) and US (SP9) standards.

Drop-test. Height in cm.



“Single-Tube” design with exclusive pressure-distribution system

The exclusive pressure-distribution system controls the airflow to the manometer during inflation to provide a smooth and accurate increase in pressure and a long manometer life. With this innovation the HEINE GAMMA G Series offers an important durability feature previously only available in “double-tube” systems, in a convenient, single-tube design.



[01] MANOMETER
[02] TUBE TO CUFF

The innovative HEINE “GAMMA G” Series latex-free sphygmomanometers sets a new standard in performance, quality and value. A functional and modern design, in combination with a high-quality precision manometer and durable construction, guarantee reliable performance over many years in any clinical environment.

GAMMA G Series / Confirmation of clinical validation:

The GAMMA G Series is clinically-validated under the European Union protocol on hypertension. The validation protocol for the GAMMA G range confirms the long-term accuracy of measurement of the systolic and diastolic blood pressure.



[HEINE GAMMA® G5]

GAMMA G5 | ORDER NO.:

WITH ADULT CUFF	M-000.09.230
WITH CHILD CUFF	M-000.09.231
PRACTICE KIT WITH CHILD,	
ADULT AND SMALL ADULT CUFF	M-000.09.555
WITH ADULT CUFF, WITHOUT	
ZIPPER POUCH, BOX OF 10	M-000.09.561

The HEINE Standard Sphygmomanometer.

The core of the GAMMA G5 is the quality and durability of the housing and spoon made of strong, impact-resistant 2-component thermoplastic. The integration of an intricate shock-absorbing system protects the manometer.

The extra large insufflation bulb is comfortable to hold, enabling fast inflation of the cuff.

- ⌘ **Latex-Free.**
- ⌘ **Shockproof.**
- ⌘ **Exclusive Pressure-Distribution System** allows swift and precise readings and ensures a long manometer life.
- ⌘ **Pressure-Proof Manometer:** Ensures durability and reliability. Corrosion-free.
- ⌘ **Micro-Filter:** Protects valve and manometer.

- ⌘ **Large Diameter Scale:** (56mm) is easy to read.
- ⌘ **Large Insufflation Bulb:** 50% larger bulb volume for fast inflation of the cuff.
- ⌘ **Precision Air-Release Valve:** Provides fine and precise air-release rate adjustment, quick deflation, and extreme durability. A rubber-stop prevents valve-lock.
- ⌘ **Housing** is moulded of shockproof plastic with shock-absorbing, anti-slip coating.
- ⌘ **Adjustable Spoon (Polyamid Thermo Plastic):** For both right- or left-handed operation.
- ⌘ **High quality, Single-Tube Cuff with Velcro Fastener.** Washable.
- ⌘ **5 Year Warranty*.**

*against defects in design, materials and workmanship for a period of 5 years from the date of purchase



[HEINE GAMMA® G7]

GAMMA G7 | ORDER NO.:

WITH ADULT CUFF	M-000.09.232
WITH CHILD CUFF	M-000.09.233
PRACTICE KIT WITH CHILD,	
ADULT AND SMALL ADULT CUFF	M-000.09.554
WITH ADULT CUFF, WITHOUT	
ZIPPER POUCH, BOX OF 10	M-000.09.560

The model with aluminium alloy construction and precision valve.

The GAMMA G7 is unique in its advanced construction, utilising a housing made of aluminium alloy and advanced thermoplastics, producing a shockproof and modern design.

Combines all design and quality features of the GAMMA G5, except:

- ⌘ **Upper Housing** made of aluminium alloy with a shock-absorbing ring protecting indicator and glass.
- ⌘ **Lower Housing** is moulded of shockproof plastic with shock-absorbing, anti-slip coating.
- ⌘ **Light Alloy, Adjustable Spoon:** For both right- or left-handed operation.
- ⌘ **10 Year Warranty*.**

*against defects in design, materials and workmanship for a period of 10 years from the date of purchase



[HEINE GAMMA® GP]

Standard Instrument with Push-Button Valve.

Equipped with an innovative, absolutely reliable and easy-to-operate push-button valve, the new GAMMA GP combines all the design and quality features of the tried and proven GAMMA G5.

- ✦ **Innovative Push-Button Valve.** Fine and precise air-release rate adjustment, ultra-fast air-release, practically wear-free. Easy to operate by both left- and right-handed users.

- ✦ **5 Year Warranty*.**

*against defects in design, materials and workmanship for a period of 5 years from the date of purchase



GAMMAGP | ORDER NO.:

WITH ADULT CUFF	M-000.09.242
WITH CHILD CUFF	M-000.09.243
PRACTICE KIT WITH CHILD,	
ADULT AND SMALL ADULT CUFF	M-000.09.556
WITH ADULT CUFF, WITHOUT	
ZIPPER POUCH, BOX OF 10	M-000.09.562

[HEINE GAMMA® XXL LF]

The versatile professional model.

The Classic GAMMA XXL LF comes in a latex-free (LF) version with a range of mounting options that makes it the ideal sphygmomanometer in any clinical setting.

The oversized, adjustable scale is easy to read, even from a distance.

Since the GAMMA XXL LF is non-mercurial, it is both safe and environmentally-friendly.

- ✦ **Latex-Free.**
- ✦ **Durable Polycarbonate Construction** will hold up in the busiest hospital and practice environments.
- ✦ **Oversized 130mm diameter scale** can be pivoted for clear, easy reading from any angle.

- ✦ **Long Spiral (3m) Inflation Tube** for increased working distance.
- ✦ **High Quality Indicator** is precisely calibrated and durable.
- ✦ **Precision Air-Release Valve** made of metal.
- ✦ **Pressure-Proof Manometer:** Ensures durability and reliability.
- ✦ **Wall, Desk, Rail and Wheeled Stand (adjustable height) Mounting Options.**
- ✦ **Mercury-free.** Environmentally friendly.
- ✦ **10 Year Warranty*.**

*against defects in design, materials and workmanship for a period of 10 years from the date of purchase

GAMMA XXL LF | ORDER NO.:

WALL MODEL LF-W	M-000.09.323
WHEELED STAND MODEL LF-S	M-000.09.326
RAIL MODEL LF-R	M-000.09.325
DESK MODEL LF-T	
(NON SPIRAL INFLATION TUBE)	M-000.09.322

[XXL LF-S]



[XXL LF-T]



[XXL LF-R]



[XXL LF-W]

[CUFFS / ACCESSORIES]

Single-tube cuffs for GAMMA G7 / G5 / GP	Infant	6 x 28	M-000.09.700
	Child	10 x 34	M-000.09.702
	Adult small	14 x 53	M-000.09.704
	Adult	14 x 58	M-000.09.619
	Adult large	17 x 70	M-000.09.706
	Thigh	20 x 86	M-000.09.615
	Cuff for GAMMA XXL LF	Infant	6 x 28
Child		10 x 34	M-000.09.703
Adult small		14 x 53	M-000.09.705
Adult		14 x 58	M-000.09.629
Adult large		17 x 70	M-000.09.631
Thigh		20 x 86	M-000.09.625
Mounting accessories and parts for GAMMA XXL LF		Cuff holder	for Desk Mount GAMMA XXL LF-T
	Desk stand	for Desk Mount GAMMA XXL LF-T	M-000.09.142
	Universal mounting plate	for rail, wall and wheeled stand GAMMA XXL LF-R / LF-W / LF-S	M-000.09.306
	Rail clamp	GAMMA XXL LF-R	M-000.09.307
	Adjustable wheeled stand	GAMMA XXL LF-S	M-000.09.145



HEINE GAMMA® STETHOSCOPES

- ❖ **Solid chest piece** and ideal weight avoid stray noise from manipulation and ensure excellent resonance. Hygienic matt chrome finish.
- ❖ **Chest piece with non-chill rims** for the comfort of the patients.
- ❖ **Ergonomic, chrome-plated binaurals.** Individually adjustable.
- ❖ **High-quality Y-piece tubing** for perfect acoustics and sound reproduction. Length 56 cm.
- ❖ **Latex-free** materials.

[HEINE GAMMA® 3.1 Pulse Stethoscope]

GAMMA 3.1 | ORDER NO.:

M-000.09.941

For fine acoustic performance.

- ❖ **Stethoscope with solid flat chest piece.** The chest piece can be easily inserted under the sphygmomanometer cuff.
- ❖ **High-quality membrane,** 44 mm dia.

- ❖ **Single-tube construction** with excellent sound transmission.
- ❖ **Includes ear olives** (large/soft) for sound insulation.
- ❖ **Optionally available:** Ear olives, (small/hard) (M-000.09.946).



[HEINE GAMMA® 3.2 Acoustic Stethoscope]

GAMMA 3.2 | ORDER NO.:

M-000.09.942

For excellent acoustic performance.

Technical data as GAMMA 3.1, but with a double chest piece.

- ❖ **Combined double chest piece** with two different frequency ranges to choose from: A high-quality membrane on one side and a conical acoustic head on the other side.

- ❖ **High-quality membrane,** 44 mm diameter for excellent acoustic performance and sound transmission.
- ❖ **Cone diameter 30 mm.**



[HEINE GAMMA® 3.3 Acoustic Stethoscope]

GAMMA 3.3 | ORDER NO.:

M-000.09.943

For the examination of children.

Technical data as GAMMA 3.2.

- ❖ **Reduced contact area** of 31 mm makes the GAMMA 3.3 ideal for examining children.



[HEINE GAMMA® C3 Cardio Stethoscope]

GAMMA C3 | ORDER NO.:


M-000.09.944

Specially-developed for cardiology.

- ❖ **Exceptional acoustic performance,** particularly in the cardiac frequency range.
- ❖ **Combined double chest piece** with two different frequency ranges. Instant selection of a double membrane on one side or a conical acoustic head on the other.

- ❖ **High-performance double membrane,** 45 mm diameter gives ideal sound transmission and perfect skin contact for exceptional sound reproduction.
- ❖ **Two-in-one tube system** for enhanced performance. Very supple, yet kink-resistant with exceptional acoustic reproduction.
- ❖ **Includes ear olives** (large/soft) and one pair of replacements (small/hard) for sound insulation.





■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES


A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet the relevant international standards (ISO/IEC) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



GERMANY

 HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7
82211 Herrsching
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com
www.heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
PO Box 7218 Warringah Mall
NSW 2100
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobeläckerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: info@heineschweiz.ch

We reserve the right to change specification without notice.
© 2016 HEINE Optotechnik GmbH & Co. KG. All rights reserved.

Your distributor

04/16. A-000.00.215 e

