



Choose the stethoscope that helps you achieve your best.

3M™ Littmann® Stethoscopes offer superb acoustic performance, consistently high quality, and an outstanding value. Nearly all Littmann stethoscopes use tunable technology, a 3M invention that gives you the ability to hear different frequency sounds without turning over and repositioning the chestpiece.

Experience the power of tunable technology.



LIGHT CONTACT:
Hold the chestpiece with light pressure to hear low frequency sounds.



FIRM PRESSURE:
Press the chestpiece a little more firmly to hear higher frequency sounds.

Lightweight		Classic				Cardiology		Electronic											
6 Acoustic Performance Ideal for taking blood pressure readings and making limited patient assessments.		7-8 Acoustic Performance Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients.				9-10 Acoustic Performance Designed to meet the needs of a demanding clinical environment. Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation.		10+ Acoustic Performance High-performance auscultation. Hear more of what you need to hear with amplification, ambient noise reduction, and optional connectivity.											
Lightweight II S.E. stethoscope		Select stethoscope		Classic II Paediatric and Infant stethoscopes		Classic II S.E. stethoscope		Classic III stethoscope		Master Classic II stethoscope		Cardiology IV stethoscope		Master Cardiology stethoscope		Electronic stethoscope Model 3100		Electronic stethoscope Model 3200	
Double sided chestpiece 118 grams Tunable diaphragm 2 year warranty Chestpiece: 5.4cm		Single sided chestpiece 120 grams Tunable diaphragm 2 year warranty Chestpiece: 4.4cm		Double sided chestpiece 105, 95 grams 3 year warranty Chestpiece: Paediatric 3.3cm, Infant 2.7cm, 1.9cm		Double sided chestpiece 135 grams Tunable diaphragm 3 year warranty Chestpiece: 4.4cm		Double sided chestpiece 150 grams Single-piece tunable diaphragm on each side of the chestpiece Convertible open bell 5 year warranty Chestpiece: Adult side 4.3cm, Paediatric side 3.3cm		Single sided chestpiece 160 grams Tunable diaphragm 3 year warranty Chestpiece: 4.4cm		Double sided chestpiece 177 grams Available in 69cm Single-piece tunable diaphragm on each side of the chestpiece Dual-lumen tubing 7 year warranty Chestpiece: Adult side 4.3cm, Pediatric side 3.3cm		Single sided chestpiece 185 grams Available in 69cm Tunable diaphragm Dual-lumen tubing 7 year warranty Chestpiece: 5.1cm		185 grams 2 year warranty Chestpiece: 5.1cm		Bluetooth® wireless technology Recording and playback 185 grams 2 year warranty Chestpiece: 5.1cm	

Bluetooth® is a registered trademark of Bluetooth SIG. 3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology IV, Master Classic II, Classic III, and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in UK. © 3M 2017. All rights reserved. 70-2011-5462-5

3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide



Model 3200 stethoscope	69 cm (27 in)	3200BK27	3200BU	3200NB																																				
Model 3100 stethoscope	69 cm (27 in)	3100BK27	3100BU	3100NB																																				
Master Cardiology stethoscope	69 cm (27 in)	2160 2161 2175 ² 2176 ⁴	2163	2164	2165	2167				2178									2182 ⁴																					
NOW AVAILABLE		6152 6162 ⁸ 6163 ⁷ 6164 ⁹ 6170 ¹¹ 6165 ¹⁰ 6177 ¹¹ 6179 ¹²	6153	6154 6168 ⁷	6155	6156 6166 ⁸				6158 6178 ⁸	6157 6169 ¹¹											6171 ⁸																	6161	
Master Classic II stethoscope	69 cm (27 in)	2144L 2141 ¹	2146	2147	2632																																			
Classic III stethoscope	69 cm (27 in)	5620 5803 ¹ 5811 ⁴ 5870 ⁵	5627 5868 ¹	5622 5867 ¹	5624	5831	5809 ³	5621	5630															5633																
Classic II S.E. stethoscope	71 cm (28 in)	2201		2205																																				
Classic II Paediatric stethoscope	71 cm (28 in)	2113		2123																																				
Classic II Infant stethoscope	71 cm (28 in)	2114																																						
Select stethoscope	71 cm (28 in)	2290	2293						2301	2294	2296	2291																												
Lightweight II S.E. stethoscope	71 cm (28 in)	2450	2451									2452		2455																										

Special chestpiece and eartube finish combinations are available for selected 3M™ Littmann® Stethoscopes. See above for tubing colour options.



1 Black Edition – black finish chestpiece and eartubes.



3 Copper Edition – copper finish chestpiece and eartubes.



5 Rainbow Edition – rainbow finish chestpiece and brass finish eartubes.



2 Brass Edition – brass finish chestpiece and eartubes.



4 Smoke Edition – smoke finish chestpiece and eartubes.



6 Rainbow Edition – rainbow finish chestpiece and black finish eartubes.

Colours shown are approximate.



3M United Kingdom PLC
Morley Street
Loughborough
Leicestershire
LE11 1EP

www.3M.co.uk/healthcare
www.littmann.co.uk

3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Cardiology IV, Master Classic II, Classic III, and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in UK. © 3M 2017. All rights reserved. 70-2011-5462-5

New colour editions for the 3M™ Littmann® Cardiology IV™ Stethoscope include special finishes on the chestpiece, eartubes and stems. See above for tubing colour options.



7 Black Edition – black finish chestpiece, stem and eartubes.



9 Brass Edition – brass finish chestpiece and smoke finish stem and eartubes.



11 Mirror Edition – mirror finish chestpiece and stem and stainless steel finish eartubes.



8 Smoke Edition – smoke finish chestpiece, stem and eartubes.



10 Rainbow Edition – rainbow finish chestpiece and black finish stem and eartubes.



12 Champagne Edition – champagne finish chestpiece and smoke finish stem and eartubes.